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OBTURATION SYSTEM SY-FILL

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Perfect Ergonomics

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Lightweight body

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Easy to operate

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Wireless operation

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Extreme temperature control

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The gutta-percha gun ensures the melting of the gutta-percha stick in a short time

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Make sure the gum doesn't get out of the root canal

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OBTURATION SYSTEM SY-FILL

Features(Obturation gun)

- Wireless 3D filling design
- Comfortable holding and easy to operate
- Accurate temperature-control and thermal protector cap can prevent accidental scald.
- Screw design can effectively prevent excessive glue.
- It can be heated to 230°C in 15s, which can melt all kinds of gutta percha.
- Large battery capacity and it can use longer time.
- Auto idle power saving mode, the device will turn off automatically after 15mins of non-use.

Features(Obturation pen)

- Wireless design, easy to operate and light body
- Adjustable hot needle size and temperature mode, gutta percha vertical condensation technique can compress gutta percha more compact.
- It has own heating protection function, the heating maximum time is 10s, more efficient and secure.
- Heating button 360° stereo operation, accurate control and 3D filling more thoroughly.
- It can be heated to 230°C in 5s.
- Optional temperature: 160/180/200/230°C

Filling technology

The filling of the apical 1/3 of the root canal using obturation pen			
Step one	Step two	Step three	Step four
Mark the working length of the pen tip, 4-7mm from the apex, to ensure the fluidity and sealin of the filling.	Insert the master gutta-percha point infiltrated with sealant.	Put the pen tip and heat it, cut it from the middle and press down to 6-9mm, stay for 1.5-2s.	Stop heating the tip with continuous pressure for 8s. Heat the tip again for 1-1.5s. Press gutta percha to 4-7mm from the heat tip.

The filling of the coronal 2/3 of the root canal using obturation gun			
Step one	Step two	Step three	Step four
Put the tip of the gun needle into root canal, in close contact with master gutta percha, heat for 3-5s	Continue to fill by squeezing the teether.	While pushing back the tip of the gun, squeeze the melted gutta percha into root canal and press the gutta percha around it.	Completely fill the entire root canal.